

Adapt Medical Adhesive Remover Spray/Adapt Universal Remover Wipe

General Information

- Helps remove adhesive skin barriers and/or tape
- Removes adhesive residue from the skin
- Silicone-based formulation
- No sting from alcohol

Using Adhesive Remover Spray

1. Shake well before using.
2. Peel back one edge of the skin barrier or tape border.
3. Hold the can so that the arrow faces the base of the skin barrier or tape, and spray until the adhesive releases. The 360° spray technology minimizes waste, and will allow you to spray from any angle.
4. Ensure your skin is dry and apply your new ostomy pouching system as usual.

Using Universal Remover Wipes

1. Tear open the package and take out the wipe.
2. Press or rub the wipe against your skin as you peel off your skin barrier and/or tape.
3. Press or rub the wipe over any areas of the skin where tape residue remains.
4. Ensure your skin is dry and apply your new ostomy pouching system as usual.



Adapt Medical Adhesive Remover Spray/Adapt Universal Remover Wipe

Stock No 7737/7760



Spray - Spray until the skin barrier releases.



Wipes - Rub the wipe over skin as you peel off skin barrier.



Wipes - Clean off remaining tape residue.

FREQUENTLY ASKED QUESTIONS

Is an adhesive remover always needed?	No, this product is optional. Use it to remove adhesive residue from your skin, or to make it easier to remove your skin barrier or tape.
Will this be uncomfortable to use?	This silicone-based product is formulated to be a no-sting alternative. When used correctly, the adhesive remover will help to make the removal of your pouching system easier and help reduce the potential for skin stripping.
Do I need to wash the adhesive remover off of my skin?	The adhesive remover evaporates quickly and dries without residue, so the next barrier can attach securely. If you wash your skin after use be sure to rinse and dry thoroughly before the next pouch application.

For product questions, sampling needs, or detailed clinical questions concerning our products in the US, call **1.888.808.7456**. In Canada call **1.800.263.7400**.

Routine follow-up with your healthcare professional is recommended.

Prior to using any ostomy products/accessories be sure to read all product inserts and labels.

The Hollister logo, Adapt and "Healthy skin. Positive outcomes." are trademarks of Hollister Incorporated.
© 2017 Hollister Incorporated. 923073-1217 US-00218

Hollister Incorporated
2000 Hollister Drive
Libertyville, Illinois 60048
1.800.323.4060

www.hollister.com

Hollister Limited
95 Mary Street
Aurora, Ontario L4G 1G3
1.800.263.7400