



VaPro Plus Pocket™

Touch Free Hydrophilic Intermittent Catheter

Protection
Discretion
Ease of Use



The new **VaPro Plus Pocket** is a hygienic, touch free hydrophilic intermittent catheter with an integrated collection bag in a smart pocket-size package.

- Catheter and collection bag are packaged in a small, discreet pouch
- The integrated 1000ml collection bag makes self-catheterisation more convenient and gives the patient more flexibility
- Patients can go about their day-to-day lives with more confidence
- Catheters are available for men, women and children
- The range is easy to use and simple to demonstrate or learn



 Hollister®

Recommend to your patients and give them back their independence

Protective tip

- Helps shield the sterile catheter, during insertion, from bacteria located within the first 15mm of the distal urethra
- Helps reduce the risk of carrying bacteria into the urinary tract

Protective sleeve

- Makes the catheter easy to grip and handle
- Provides a barrier to the catheter to help reduce bacterial contamination

Collection bag

- 1000ml Integrated collection bag for convenient catheterisation at home or on the go
- Anti-reflux valve reduces risk of urine backflow

Designed for daily life

- Full-length (20cm Female, 40cm Male) catheter
- Pure catheter hydration technology makes it ready for use right out of the package
- Designed to be evenly lubricated for easier insertion and withdrawal
- Easy-open hole and tear system is ideal for users with limited dexterity

Women

Female Catheter

Pocket Size:

15.5cm x 9.6cm

Catheter length:

20cm



Men

Male Catheter

Pocket Size:

15.5cm x 10.6cm

Catheter length:

40cm



Stock No.	Length (inches)	Size (Charriere)	Box Qty
71082	20cm (female)	Ch 8	30
71102	20cm (female)	Ch 10	30
71122	20cm (female)	Ch 12	30
71142	20cm (female)	Ch 14	30
71084	40cm (male)	Ch 8	30
71104	40cm (male)	Ch 10	30
71124	40cm (male)	Ch 12	30
71144	40cm (male)	Ch 14	30
71164	40cm (male)	Ch 16	30

This intermittent catheter is a flexible, tubular device that is inserted through the urethra to drain urine from the bladder. **WARNING:** To help reduce the potential for infection and/or other complications, do not reuse. If discomfort or any sign of trauma occurs, discontinue use immediately and consult your healthcare practitioner. **PRECAUTION:** Please consult your healthcare practitioner before using this product if any of the following conditions are present: • Severed urethra • Unexplained urethral bleeding • Pronounced stricture • False Passage • Urethritis – inflammation of the urethra • Prostatitis – inflammation of the prostate gland. Self-catheterisation should follow the plan of care and advice given by your healthcare practitioner and be carried out only in accordance with the instructions for use provided. Because catheterisation frequency varies by person, the recommended frequency of your catheterisation should be provided by your healthcare practitioner. For any other questions about your catheterisation, contact your healthcare practitioner.

Please visit hollister.co.uk or call 0800 521 377 to request a sample



Hollister Limited – UK

Rectory Court, 42 Broad Street
Wokingham, Berkshire RG40 1AB